



**LC-38**

**2.0MP HDMI Digital Camera**



## LC-38 2.0MP HDMI Digital Camera

LC-38 2.0MP HDMI Digital Camera are economic version with simple and compact structure CMOS eyepiece cameras. USB Type C is used as the data transfer interface. Also the LC-38 2.0MP HDMI Digital Camera comes with advanced video & image processing application KoPa Capture.

### FEATURES

- Ultra-wide dynamic and ultra-low-light camera.
- HDMI 1080P 60Hz preview, 30FPS record video.
- USB high-speed video transmission.
- Rich PC software functions.
- Built-in multifunctional image processing software.
- Local storage supports up to 128G and supports hot plug.
- Equipped with a wired controller as standard accessory to reduce work intensity.

### SPECIFICATION

Model	LC-38
Physical resolution	2.0MP
Sensor model	SONY CMOS IMX291
Exposure mode	Roller shutter exposure
Highest resolution	1920 x 1080 (2,073,600)
Sensor size	1/2.8"
Pixel size	2.9 $\mu\text{m}$ x 2.9 $\mu\text{m}$
Dynamic Range	128dB
SNR	30dB
Spectral response	380~650nm
Exposure	Real-time auto, single auto, manual adjustment
White balance	Real-time automatic, single automatic, manual R G B adjustment
Preview resolution	HDMI: Auto-adaptive display, up to 1920X1080i60, 1020X1080P30
Record format	"Snapshot Format: JPG Resolution: HDMI 1920x1080, Record Format: MOV Resolution: HDMI 1920X1080P30,1280x720P60"
Record path	TF card (up to 128G,hot pluggable)
Power supply	USB TYPE C 5.0V DC (maximum current 1A)
Image output interface	HDMI TYPE C(Mini HDMI)

DIAGRAM

